



MILFORD QUADBOX

The solution for storage and transport of liquids

About Milford IBC

Milford has been established for more than forty years, and during this period has supplied an increasing customer base comprising local companies, as well as **national and international organisations**.

Milford continues to supply these companies today, and welcomes the opportunity to **create new equipment** solutions. These solutions are developed with the needs of the client in mind, and for this reason solutions differ vastly from project to project. Common outcomes are gains in processing efficiency, and **optimised transport costs**.

These logistics equipment solutions can be applied across many industries including:

- Post and logistics
- Food manufacturing and processing
- Beverages, juice and concentrates
- Dairy products
- Automotive logistics
- Chemicals and industrial products

Reusable folding IBC containers are 100% recyclable, and **reduce the need for single use packaging** which is good news for the environment. They are designed and built to deliver many years of trouble free operation with minimal need for servicing.

When these IBCs are folded they occupy significantly **less space in your warehouse**, and deliver significant **savings in return transport**.

Milford QuadBox details

The Milford QuadBox has been designed as a cost efficient reusable folding IBC to be used with commonly available single use liner bags or bladders. It can be used for the storage and transportation of a variety of liquids.

The side panels are formed from tubular steel frames then fitted with long lasting UV stabilised food grade polypropylene plastic panels. The Milford QuadBox is finished in a durable anti corrosion finish to ensure many years of trouble free operation.

Assembly and folding the Milford QuadBox is quick and easy, and ticks all the boxes for ergonomics. The Milford QuadBox also features a graduated window to allow checking on the contents during filling or emptying, and a low level liquid outlet to accept all common 2" discharge and filling glands. This outlet can be sealed to ensure tamper evidence during transport and storage.

Milford QuadBox specifications

- Base frame size is 1140mm x 1140mm
- Height is 1185mm assembled and 385mm folded
- Internal capacity is 1,100 litres or 1,500 kg
- The Milford QuadBox can be stacked
 - Assembled - 2 high in transport or 5 high static
 - Folded - 6 high in transport or 14 high static
- Four way entry strengthened pallet base
- Easy disassembly to ensure thorough cleaning
- Low level discharge for efficient emptying

Milford IBC UK Pty Ltd
Gloucester House
399 Silbury Boulevard,
Milton Keynes, MK9 2HL

T: +44 (0) 1908 303 700
F: +44 (0) 1908 678 965
E: salesuk@milfordibc.com
www.milfordibc.com

